

Manual Operation Theatre Table For Operating Room

Introduction


1.3001 is hydraulic mechanical operating table. The table up-down is hydraulic controlled.

2.The column and base cover are all stainless steel 304#, which is convenient for clean and disinfection.



Main Specification

Main Specification	
Length/Width	2000mm/480mm
Height	720mm-980mm
Trendelenburg/ Revers Trendelenburg	$\geq 30^\circ / \geq 20^\circ$
Lateral Tilt	$\geq 18^\circ / \geq 18^\circ$
Head Section Up/Down	$\geq 45^\circ / \geq 90^\circ$
Leg Section Up/ Down/ Outside	$\geq 15^\circ / \geq 90^\circ / \geq 90^\circ$
Back Section Up/ Down	$\geq 75^\circ / \geq 18^\circ$
Kidney Bridge Lift	115mm
Load Capacity	250kg



Tel. +90 850 228 87 00

info@biostar-medikal.com.tr

www.biostar-medikal.com.tr

