

3001S Manual Operation Theatre Table For Operating Room

Introduction

This operating table is suitable for surgical procedures of the head, neck, thorax, abdomen perineum and extremities as well as for gynecological, obstetrical, urological, otorhino laryngological and cranial operations. It can be raised, lowered, Reversed Trendelenburg, Trendelenburg, Lateral Tilt. Leg section can be rotated, opened, disassembled and adjusted easily. It is most adapted to the urological operation. The table top made of high strength plastic can move backwards and forwards. **X-ray examination can be carried out on it.**

It can be raised or lowered hydraulically by means of an oil pump. Various positions can be controlled from both side. Chassis inside outside cover of this table are stainless steel.



Main Specification

Main Specification	
Length/Width	2000mm/480mm
Height	720mm-980mm
Trendelenburg/ Revers Trendelenburg	$\geq 30^\circ / \geq 20^\circ$
Lateral Tilt	$\geq 18^\circ / \geq 18^\circ$
Head Section Up/Down	$\geq 45^\circ / \geq 90^\circ$
Leg Section Up/ Down/ Outside	$\geq 15^\circ \geq 90^\circ / \geq 90^\circ$
Back Section Up/ Down	$\geq 75^\circ / \geq 18^\circ$
Kidney Bridge Lift	115mm
Load Capacity	250kg



Tel. +90 850 228 87 00

info@biostar-medikal.com.tr www.biostar-medikal.com.tr

