F3 **Fetal Monitor**

Linn



About Edan

Edan is a healthcare company dedicated to improving the human condition around the world by delivering value-driven, innovative and high-quality medical products and services. For over 20 years, Edan has been pioneering a comprehensive line of medical solutions that address a broad range of healthcare practices including:

- Diagnostic ECG
- Ultrasound Imaging
- In-Vitro Diagnostics

- Patient Monitoring • OB/GYN
- Point-of-Care Testing
- Veterinary

Healthcare professionals around the world depend on Edan's breakthrough medical technologies and outstanding customer support.



Global Headquarters: Edan Instruments, Inc. | 15 Jinhui Road, Pingshan District, Shenzhen, 518122 P.R. China | +86.755.26898326 | www.edan.com | info@edan.com U.S. and Canada inquiries: EDAN Diagnostics, Inc. | 9918 Via Pasar, San Diego, CA 92126 +1.858.750.3066 www.edandiagnostics.com edan-info@edandiagnostics.com

© Edan Instruments, Inc. All rights reserved. Features and specifications are subject to change without prior notice. No reproduction, copy or transmission may be made without written permission. Not all products or features are available in all countries, contact Edan for local availability.

CE0123

ENG-OB & GYN-F3-V3.5-20200101









Multi-display Modes Multi-display Modes Switching Switching</l

Functional keys and knob provide various shortcuts to achieve functions for clinical use. The 'start' button could be configured to integrate patient info and printing. This could help the doctor to simplify the workflow and work with 1 button only.

USB Storage Capacity



F3 supports USB interface, providing house call service.

> Data Management





Based on our own design language, F3 gets an iF design award for its compact, lightweight and ergonomic design.

With its elegant design, bright color screen, and enhanced functionality, F3 fetal monitor would address the needs of obstertrical departments in both clinics and hospitals. F3 is amazingly lightweight and portable, which makes it suitable for outpatient service. With the latest technologies from EDAN, F3 offers an extensive set of monitoring parameters such as FHR, TOCO, DECG, IUP and fetal movement.

Convenient to Carry or Transport

- Compact, portable and lightweight
- Long life Li-ion battery for outpatient service

Easy to Read & Operate

- 5.6" foldable color TFT screen
- Large numeric and waveform display
- Simple workflow with silicone buttons for operations

Advanced FHR Detection Technology

- Signals Overlap Verification to differentiate twins FHR
- FHR signal quality indicator helps optimize the probe position

Switching among three display modes, user may choose the most suitable one for clinical use.



F3 supports USB interface, providing enlarged storage capacity, which is a quite good solution for

With the help of EDAN Insight software, doctors can transmit CTG data from F3 to PC for data management.

Flexible Printing Options

- Built-in thermo printer
- Print on 150/152 mm width recording paper
- 1/2/3 cm/min real-time printing speed
- 15 mm/sec fast printing speed for history traces

Powerful Data Management

- 60 hours built-up memory for seamless monitoring
- MFM-CNS network for remote monitoring
- MFM-CNS Lite for data management on PC
- USB port to enlarge storage capacity