

## Electric Suction Apparatus



**Features:** Adopting oil-free piston pump

Clean-without pollution of oil and smoke

Convent-no need add oil to maintain daily

Safe-no plus pressure during using

No flow backwards pressure when the machine stops,so the liquid won't flow backwards

Complete plastic panels, hand-switch and foot-switch for operating easily. Select negative pressure at will with no level for pressure adjustment, and overflow protection mechanism, safe and reliable.

### Principle technical data:

**Max negative pressure:**  $\geq 0.09\text{MPa}$

**Noise:**  $\leq 65\text{ dB (A)}$

**Power:**  $\leq 180\text{VA}$

**Pumping rate:**  $\geq 20\text{L/min}$

**Reservoir capacity:** 2500mL/pc, 2pieces

**Note:** Not suitable for continuous use in the induced abortions and for the inflammable or explodable surroundings. Max continuous working time is 30 minutes, continuity rate: 50% Electrical appliance safety requirement: Class I, Type B

